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
You have been prescribed a FroggyMouth

This device corrects swallow function by creating a new pattern: during the swallow, the teeth come lightly in contact and the tongue pushes against the palate.

Your provider has given you a FroggyMouth


Liberate yourself of bad habits

With a Very Simple Approach

In Just 15 minutes a Day while watching TV

All of these clinical conditions have one thing in common: they are all related, in different ways, to tongue posture. In fact, for 40% of the population, what is known as a tongue thrust may still be present even after the first molars come in. This leads to poor tongue function and disrupted orofacial functions that eventually impact optimal growth.

By correcting swallow dysfunction, FroggyMouth helps by:

- 🕒 Shortening the **duration of the treatment**
- 🏠 Improving the **quality of results**
- 🌟 Reducing the **risk of relapse**

How does FroggyMouth work ?

The FroggyMouth should be worn 15 minutes a day, in front of a TV. While the child wears the device, the swallow pattern changes automatically without disrupting sleep quality, which is very important for children.

It is mandatory to follow accurately **FroggyMouth Rules of Use**

Wear it for 15 minutes a day to quickly and effectively develop a new way of swallowing. Routine is essential.

Wear in front of a raised TV screen or computer. This ensures that the horizontal visual plane is intact thus preventing from any posture with a bowed head.

The patient's attention should be held by the screen so that the new way of swallowing becomes automatic and natural.

How to wear FroggyMouth ?